


## SP-2 Capsule Separating Machine



### Instructions:

SP-2 Capsule Separating Machine(Decapsulator) is new equipment for emptying all kinds of abnormal capsule products and recovering the materials quickly,so that the materials can be reused in capsule filling machine.

### Competitive Advantages:

- 1.capsule cap is completely separated, the capsule no deformation and capsule fragments are avoided.
- 2.the machine can be applied to the capsule of all models and all the materials without changing molds, nor adjusting equipments.
- 3.the machine can be applied to abnormal capsules of various shapes, such as deformation, concaved cap,fusion,adhesion and so on.
- 4.the machine can be applied to contents of all forms,such as powder, pellets, granules,pills,etc,without imposing any impact properties of contents.
- 5.the work efficiency can reach 4000 to 12000 units per minute.
- 6.for capsules of all types and unusual forms, the machine can reach the retuning rate of over 99.9%.
- 7.the machine performed a very low failure rate, and can be operated flexibly and conveniently. it is easy for cleaning and maintenance.

SP-2 Decapsulator is used in more and more pharmaceutical companies,it works together with capsule fill machine, the cooperative operation of both can further improve the yield of capsule



---

products and provides more effective guarantee for product quality. it also plays an important role in saving costs for customers .SP-2 is regarded as "The best partner of capsule filling machine"by customers.

### **Technical Parameters:**

Model	SP-2
Capsule Size	00#---5# Capsule
Production Capacity	Max. 216 000 pcs/h
Air Consumption	5-7 Kg f/cm <sup>2</sup>
Power	220V/50Hz/1.5KW
Packing Size	410×510×860 mm
Packing Weight	130 Kg